


Mini-Mold™

Light cure your own aesthetic attachments in seconds!

Simple 2-Step Technique



STEP 1

Fill the selected Mini-Mold Tip with light cure resin.



STEP 2

Place the Mini-Mold Tip on desired tooth and light cure.



- Invisible aesthetic attachments
- As fast as 10-second cure time
- Saves chair time and money
- Safe for sensitive patients

Just fill the tip, place on tooth, and cure. In as little as 10 seconds the Mini-Mold Tip precisely casts a light-cured resin attachment.

Choose From 7 Different Mini-Mold Treatment Aids

Each Mini-Mold Tips Treatment Aid Includes:

10 Mini-Mold Tips (same type), 1 tip handle, and detailed step-by-step instructions (Bite-Ramp 8 per pack)



Bite-Ramp Mini-Mold Tips

Creates a bite-opening ramp on the lingual of the upper anterior teeth.

Product Code: 51392

Price £37.50 pack



Bracket Mini-Mold Tips

Creates edgewise brackets, ideal for minor tooth adjustments.

Product Code: 51396

Price £37.50 pack



Button Mini-Mold Tips

Creates a functional low profile button quickly and easily.

Product Code: 51394

Price £37.50 pack



Habit-Control Mini-Mold Tips

The small, conical spike, ends thumb sucking and tongue thrusting.

Product Code: 51395

Price £37.50 pack



Wire Bonder Mini-Mold Tips

Ideal for making direct bond lingual retainers.

Large - 4.5mm Diameter

Product Code: 51393

Small - 3mm Diameter

Product Code: 51397

for wire sizes up to .020"

Price £37.50 pack




Tube Bonder Mini-Mold Tips

Designed to create a temporary tube.

Product Code: 51391

Price £37.50 pack



Posterior Molar Pad Mini-Mold Tips

Designed to open the bite from the posterior teeth.

Product Code: 51390

Price £37.50 pack



Mini-Mold Starter Kit

Product Code: 51399

Kit includes: 24 Mini-Mold Tips. 3 of each style, 2 tip handles, and detailed step-by-step instructions

Price£85.00 kit

First Time Order? Try our Starter Kit, it's a great way to sample every style of Mini-Mold.