

Quick-Fit™

Ortho Technology Molar Bands



Micro-Etched Interior Surface
Ensures maximum retention

Rolled Edges
No sharp edges
maximum patient comfort

Permanent Laser Marking
For easy identification

Anatomically Contoured
Fast placement and accurate fit

Straight Interproximal Wall
Makes band seating easy

Ideally Tempered
For easy contouring

Many more band options are available

Quick-Fit Bands

Prewelded Molar Bands

- Price £3.50 each
- 100 bands £3.25 each
- 200 bands £3.15 each



Maxillary Combination Convertible

Occlusal Headgear Tube
Assembly with Lingual Cleat

Standard Preweld Position

Buccal Attachments to Molar Bands

The main archwire slot is centred occlusal/gingivally, and the hook is always placed towards the gingiva. Standard angulation is always 0 degrees, which means the main slot runs parallel with the occlusal margin of the band. The mesial face of the tube bisects the buccal cusp on the maxillary arch, while on the mandibular arch the main body of the tube is centred on the buccal indent.

Rx	Slot Torque	Distal Offset	Headgear Tube	Slot Size	Left	Right
Roth	-10°	8°	.045"	.022"	VC6-34-206-000	VC6-34-205-000
MBT	-14°	10°	.045"	.022"	VC6-34-018-000	VC6-34-017-000

Lingual Attachments to Molar Bands

Standard positioning for lingual attachments are as follows for each attachment:
Lingual Sheaths, Cleats, and Lugs are centred mesial/distal as well as occlusal/gingival.
Buttons are centred on the mesial cusp and are offset 1mm towards the occlusal margin from the centre of the band.



Mandibular Single Convertible

Assembly with Lingual Cleat

Rx	Slot Torque	Distal Offset	Slot Size	Left	Right
Roth	-25°	5°	.022"	VC6-36-208-000	VC6-36-207-000
MBT	-20°	0°	.022"	VC6-34-062-000	VC6-34-061-000