


Engage-R Tool



The dual tips on one end span the In-Ovation R bracket to grab the arch wire, allowing you to easily snap the clip closed. The other end has a small ball that fits perfectly into the sliding cap for easy opening.

Product Code: ODG88005

Price £49.95 each

Engage-C Tool



The dual tips on one end span the In-Ovation C bracket to grab the arch wire, allowing you to easily snap the clip closed. The other end has a small "cup" that fits perfectly over the ball of the sliding clip for easy opening.

Product Code: ODG88008

Price £49.95 each

Engage-L Tool



The fork end helps push the wire into the slot of the In-Ovation L bracket. The other end is used to rotate the bracket clip open with ease.

Product Code: ODG88009

Price £49.95 each

Beaver Tail



The Beaver Tail has a flat end for opening the In-Ovation R clip and a forked end for seating the wire before closing the clip.

Product Code: ODG88001

Price £39.95 each

TN3 Bracket Positioner




The TN3 has on one end a "hockey stick" design with laser etched scribe lines for precise placement of brackets; and on the other end, a scaler for making clean-up of flash easy.

Product Code: ODG520

Price £23.85 each

Lotus SLB Explorer



To open the Lotus bracket clip, it is best to use the Explorer Hand Instrument.

Product Code: OT-2106

Heat sterilisable up to 370°F/188°C

Price £11.95 each

5 or more £9.95 each

Lotus SLB Bracket Tweezers



To close the Lotus bracket clip use the Lotus Bracket Tweezers.

Product Code: OT-202

Heat sterilisable up to 370°F/188°C

Price £14.95 each

5 or more £10.95 each