

Interproximal Stripping

Safe and Precise Interproximal Reduction

Intensiv Ortho Strip System

The Intensiv Ortho-Strip System combines a Kavo handpiece with Intensiv Ortho-Strips to provide the safest and most precise option available for interproximal reduction. It offers outstanding reliability. The Ortho-Strip kit is an irrigated system reducing friction and heat. Contouring is achieved with the flexible diamond strips and the reduction area may be approached from multiple angles as the handpiece head rotates 360 degrees. The kit includes a Kavo Contra Angle Handpiece with Head, and both double and single sided Ortho-Strips.

- Rapid and controlled reduction of enamel
- Safe and efficient
- Outstanding reliability
- Fully autoclavable strips

See page 131 for Tooth Mousse to help restore lost minerals and protect from decay & decalcification



Intensiv-Ortho Stripping Kit
Product Code: KIT97-990-00
Price£998.50 each

Code	Strips - Sold individually	Price
97-990-19	OS Double 15µm Yellow	£19.95
97-990-18	OS Double 25µm White	£19.95
97-990-17	OS Double 40µm Red	£19.95
97-990-06	OS Double 60µm Brown	£19.95
97-990-16	OS Double 90µm Clear	£19.95
97-990-12	Double Sided Set (Hack saw) (OS60/OS40/OS25/OS15)	£78.75
97-990-13	Single Sided Set (Cake Knife) (Proxoshape set, 10 blades)	£156.45
97-990-20	Proxo Planer Complete Set (Manual handpiece & 6 files - set)	£126.00
97-990-15	Proxo Planer File Set (6 files, 1 of ea. DS & SS med/fine/x-fine)	£99.75



“ I have used the Intensiv Ortho-Strip Kit for nearly [4] years and it has become one of my favourite instruments, particularly in the care of adult patients most of whom have some triangular black spaces. When used in conjunction with a fine disc or diamond to free the contact point, plus a set of measuring gauges, it gives tremendous precision in sculpting contact points. The adjustments required as part of the Invisalign plans have become very quick, easy and comfortable for the patient. I hardly pass a day without using it. ”

Dr Peter Huntley, Quality Orthodontics, Solihull



Intensiv ProxoStrip

A unique, diamond-coated metal strip for final polishing of proximal surfaces

Practical, fast, ergonomic!

A metal strip with grip secure strip ends, diamond-coated on one side, with two working areas of different grit sizes.

- Sterilisable, reusable
- Finishing and polishing in one process
- Ergonomic grip for simple, safe use
- Short strip for improved finger support while stripping



Pack of 6TOC-PX4015-6

Price£42.65 pack

Pack of 12 ...TOC-PX4015-12

Price£68.50 pack



Ideal Interproximal Reduction System

The Ideal Interproximal System has a W&H handpiece combined with the Orthofile System for safe and efficient stripping. The ultra-thin metal strip and plastic handle allow the strips to be flexible for easy contouring and they are colour coded for easy identification. The Ideal strips are sold in packs of 7.

- Safe reduction
- Flexible strips, sterilisable
- Colour coding for easy identification

Ideal Interproximal Reduction Kit

Product Code: KIT04-100-10

Price£595.00 each

Code	Strips - Sold in packs of 7	Price
04-200-15	15µm White (xx-fine)	£112.00
04-200-25	25µm Yellow (x-fine)	£112.00
04-200-40	40µm Red (fine)	£112.00
04-200-60	60µm Grey (medium)	£112.00
04-200-90	90µm Green (coarse)	£112.00
04-200-120	120µm Black (x-coarse)	£112.00
04-200-150	150µm Blue (xx-coarse)	£112.00
04-200-10	Ideal Strips pk 7 one of each of the above	£112.00

See page 131 for Tooth Mousse to help restore lost minerals and protect from decay & decalcification



SWISS MADE

90µm Green

Used for maximum reduction at the initiation of interproximal stripping and contouring. Will reduce restorative enamel, restorative materials and composites.

60µm Grey

Used for initial reduction when there is a degree of crowding or when there is existing space. Will reduce and contour enamel, restorative materials and composites.

40µm Red

Used for reduction when crowding exists or where you would like to contour to a fairly smooth finish. Will reduce enamel, restorative materials and composites.

25µm Yellow

Used for initial stripping in crowded cases and where you are aiming for a prepolished surface finish.

15µm White

Used for initial stripping in severely crowded cases and where you are finishing enamel, composites or ceramics.

Now available

120µm Black - Used for maximum reduction at the initiation of the interproximal stripping.

150µm Blue - Used for initial reduction when a larger reduction is indicated.

Gauges


Incremental Thickness Gauge

- Six stainless steel gauges
- Dimensions are 0.1mm, 0.2mm, 0.25mm, 0.3mm, 0.4mm, 0.5mm



Product Code: 25-990-25

Price£47.00 each



ARS Space Measuring Gauge

- Indicates precise measurement of ARS sites
- Also doubles as a periodontal pocket measuring instrument
- Laser markings at 3, 5 and 7mm at the tip of the instrument
- Also excellent for tucking ligature wires

Product Code: 25-500-02

Price£39.95 each